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
1. Linkow L†, Charchee R†. Theories and Techniques of Oral Implantology. St. Louis, MO: Mosby. 1970; pp. 74-76.
2. Brånemark PI†, Hansson BO†, Adell R†, et al. Osseointegrated implants in the treatment of the edentulous jaw: experience from a 10-year period. *Scand J Plast Reconstr Surg.* 1977;2:1-132.
3. Albrektsson T, Zarb GA†, Worthington P, et al. The long-term efficacy of currently used dental implants: a review and proposed criteria for success. *Int J Oral Maxillofac Implants.* 1986;1:11-25.
4. Schnitman PA†, Wohrle PS†, Rubenstein JE†. Immediate fixed interim prostheses supported by two-stage threaded implants; methodology and results. *J Oral Implantol.* 1990;16:96-105.
5. Schnitman PA†, Wohrle PS†, Rubenstein JE†, et al. Ten-year results for Brånemark implants immediately loaded with fixed prostheses at implant placement. *Int J Oral Maxillofac Implants.* 1997;12:495-503.
6. Tarnow D†, Emtiaz S†, Classi A†. Immediate loading of threaded implant at stage I surgery in edentulous arches: ten consecutive case reports with 1- to 5-year data. *Int J Oral Maxillofac Implants.* 1997;12:319-324.
7. Cooper LF†, Rahman A†, Moriarty J†, et al. Immediate mandibular rehabilitation with endosseous implants: simultaneous extraction, implant placement and loading. *Int J Oral Maxillofac Implants.* 2002;17:517-525.
8. Horiuchi K†, Uchida H†, Yamamoto K†. Immediate loading of Brånemark system implants following placement in edentulous patients: a clinical report. *Int J Oral Maxillofac Implants.* 2000;15:824-830.
9. Grunder U†. Immediate functional loading of immediate implants in edentulous arches: two-year results. *Int J Perio Rest Dent.* 2001;6:545-551.
10. Testori T†, Meltzer A†, Del Fabbro M†, et al. Immediate occlusal loading of OSSEOTITE® Implants in the lower edentulous jaw. A multicenter prospective study. *Clin Oral Impl Res.* 2004;15:278-284.
11. Pieri F†, Aldini NN†, Fini M†, Corinaldesi G†. Immediate occlusal loading of immediately placed implants supporting fixed restorations in completely edentulous arches: a 1-year prospective pilot study. *J Periodontol.* 2009;80(3):411-21.
12. Jaffin RA†, Berman CL†. The excessive loss of Brånemark fixtures in type IV bone: a 5-years analysis. *J Periodontol.* 1991;62:2-4.
13. Esposito M†, Hirsch JM†, Lekholm U†, et al. Biological factors contributing to failures of osseointegrated oral implants. (II) Etiopathogenesis. 1998;106:721-764.
14. Lekholm U†, Zarb GA†. Patient selection and preparation. In Brånemark PI†, Zarb GA†, Albrektsson T (eds): *Tissue-Integrated Prostheses: Osseointegration in clinical dentistry.* Chicago, Quintessence. 1985; pp. 199-202.
15. Friberg B†, Jemt T†, Lekholm U†. Early failures in 4,641 consecutively placed Brånemark dental implants. A study from stage I surgery to the connection of completed prostheses. *Int J Oral Maxillofac Implants.* 1991;6:142-146.
16. Adrianssens P†, Herman M†. Immediate implant function in the anterior maxilla: a surgical technique to enhance primary stability for Brånemark MKIII and MKIV implants. A randomized, prospective clinical study at the 1-year follow-up. *Appl Osseoint Res.* 2001;2:17-21.
17. Van Den Bogaerde L†, Pedretti G†, Dellacasa P†, et al. Early function of splinted implants in maxillas and posterior mandibles using Brånemark system machined-surface implants: an 18-month prospective clinical multicenter study. *Clin Implant Dent Relat Res.* 2003;5(Suppl 1):21-28.
18. Del Fabbro M†, Testori T†, Francetti L†, et al. Systematic review of survival rates for immediately loaded dental implants. *Int J Periodontics Restorative Dent.* 2006;26:249-263.
19. Calandriello R†, Tomatis M†. Simplified treatment of the atrophic posterior maxilla via immediate/early function and tilted implants: a prospective 1-year clinical study. *Clin Implant Dent Rel Res.* 2005;7:1-12.
20. Malo P†, Rangert B†, Nobre M†. All-on-4 immediate function concept with Brånemark system implants for completely edentulous maxillae: a 1-year retrospective clinical study. *Clin Implant Dent Relat Res.* 2005;7 (Suppl 1):S88-S94.
21. Balshi SF†, Wolfinger G†, Balshi TJ†. A prospective study immediate functional loading, following the Teeth in a Day protocol: a case series of 55 consecutive edentulous maxillas. *Clin Implant Dent Relat Res.* 2005;7:24-31.
22. Capelli M†, Zuffetti F†, Del Fabbro M†, Testori T†. Immediate rehabilitation of the completely edentulous jaw with fixed prostheses supported by either upright or tilted implants: a multicenter clinical study. *Int J Oral Maxillofac Implants.* 2007;22:639-644.
23. Testori T†, Del Fabbro M†, Capelli M†, Zuffetti F†, Francetti L†, Weinstein RL†. Immediate occlusal loading and tilted implants for the rehabilitation of the atrophic edentulous maxilla: 1-year interim results of a multicenter prospective study. *Clin Oral Implants Res.* 2008;19:227-232.
24. Fazi G†, Tellini S†, Vangi D†, Branchi R†. Three-dimensional finite element analysis of different implant configurations for a mandibular fixed prosthesis. *Int J Oral Maxillofac Implants.* 2011;26:752-759.
25. Baghai Naini R†, Nokar S†, Borghei H†, Alikhasi M†. Tilted or parallel implant placement in the completely edentulous mandible? A three-dimensional finite element analysis. *Int J Oral Maxillofac Implants.* 2011;26:776-781.
26. Romanos GE†, Nentwig GH†. Immediate functional loading in the maxilla using implants with platform switching: five-year results. *Int J Oral Maxillofac Implants.* 2009;24(6):1106-12.
27. English C. E. The critical A-P Spread. *The Implant Society.* March-April 1990;1(1):2-3.
28. Drago C†, Lazzara R†. Guidelines for implant abutment selection for partially edentulous patients. *Comp Cont Dent Ed.* 2010;31:14-28.

†The indicated clinicians have financial relationships with BIOMET 3i LLC resulting from speaking engagements, consulting engagements and other retained services.

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